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## Recommendations on the Automated Cleaning and Disinfection of Rigid Instruments for Minimally Invasive Surgery

*Study Group Hygiene in Minimally Invasive Surgery<sup>1</sup>*

*Our study group is composed of manufacturers of instruments, cleaning and disinfection agents, automatic washer/disinfectors and decontaminators, as well as hospital infection control experts, Central Service staff, and surgeons. Since these surgeons specialise mainly in the areas of laparoscopy and thoracoscopy, the following recommendations are formulated mainly for the instruments used in those specialties. It is especially in these procedures, as well as in gynaecology, that the instruments are contaminated internally as well as externally owing to insufflation pressure.*

*These recommendations are essentially applicable to the instruments used in other disciplines as well (otorhinolaryngology, neurosurgery, arthroscopy, etc.), since the instruments employed there are of a similar design. To be sure, the instruments are contaminated differently in those fields because of the lack of insufflation pressure (e.g., neurosurgery) or because the instruments are flushed with an irrigation fluid during use (e.g., arthroscopy).*

*According to the Medical Devices Act, surgical instruments must be prepared with a validated procedure. This requirement can be met only by automated reprocessing where the flow of the cleaning solution is assured and verified. Manual processing cannot be validated.*

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### Scope

A distinction is made between optical systems and instrument proper.

The optical system consists of rigid lenses, universal light cord, and the camera.

Instruments comprise tubing systems, cables, trocars, forceps, scissors, etc., and their attachments.

The recommendations presented here apply to the instrument proper only.

### Classification of Instrument Categories

1. Disposable instruments cannot be reprocessed.
2. Instruments which cannot be disassembled and which do not feature an irrigation adapter may no longer be used, since effective cleaning and sterilisation of their interior cannot be assured. They must be treated as disposable instruments.
3. Instruments with an irrigation adapter (female Luer-lock recommended) offer the possibility of internal cleaning. The cleaning process must be validated. The effective flow through the instrument must be assured. The hydromechanical properties of the instrument must not impede cleaning. Dead spaces must be avoided. It is desirable for the jet of cleaning solution to impinge on the joint area.
4. Take-apart instruments must be disassembled for cleaning. "Take-apart" means that they can be disassembled into the separate parts according to the manufacturer's instructions. The disassembly must be designed so as to improve the cleaning. The cleaning process must be validated and the effective flow through the instrument must be assured.

### Instrument Requirements

The composition of the materials used in the instruments, including the lenses, should allow for steam sterilisation. Reference is made to the applicable DIN standards (specifications) with respect to materials. Instruments should be designed such that all surfaces, both interior and exterior, can come into contact with the cleaning or disinfecting solution assuring its effectiveness. The surfaces must allow the instruments to be cleaned well. New instruments from the factory must be shipped in a manner ensuring that they can be prepared by standard reprocessing procedures. It must be possible to subject the instruments to steam sterilisation in their assembled state.

### Instrument Care and Cleaning in the Operating Room

In order to make processing easier, the following measures should be taken at the instrument table during and immediately after the operation:

1. Instruments that have an irrigation channel should be flushed by the instrument personnel during a break in the operation when the instrument is not

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in use (e.g. with a syringe and sterile rinsing solution). Spattering can be avoided by placing the tip of the instrument in a bottle or by immersing it in a rinsing solution. At any rate, this pre-cleaning must be effected immediately after the operation.

2. The coagulation surfaces of high-frequency instruments should be wiped off with a pad between uses during the operation. Encrustations can be dissolved with 3% H<sub>2</sub>O<sub>2</sub> solution. Before reuse, the surfaces have to be wiped again. In any case, this preliminary cleaning must be done immediately after the operation.
3. Gross soil adhering to the instruments, such as residual tissue, must be removed immediately after the operation. Reusable tubing systems must be flushed and disassembled. Instruments must be disassembled immediately after the operation, if possible by the instrument personnel, and mounted on a suitable fixture. (Never simply cast off instruments because they could be seriously damaged). The instruments must be subjected to further processing without delay. The fixture must protect the disassembled instruments from damage during transportation. If this is not guaranteed, the instrument can also be disassembled in the Central Sterile Supply Department. Immersion during transportation is principally rejected.

### Processing

MIS instruments must not be cleaned in the surgical unit. (The surgical unit comprises the operating room, anaesthesia room and recovery room, and the wash-room. A surgical department comprises one or more surgical units, as well as storage rooms and access rooms.)

Handling of used material should be limited to the necessary disassembly and connection to the automatic washer/disinfectors. Patency of instruments of Category 3 must be checked manually or automatically. Reprocessing immediately after surgery improves the cleaning results for the MIS instruments. These instruments have priority over other instruments for processing. Preparation within 8 hours is strongly recommended.

### Processing Automats

The instruments must be connected to the ports provided for them in the washer/disinfectors. The four effective components in processing, mechanics, temperature, chemistry, and time, are interdependent and must be adjusted accordingly. The processing procedure must be validated, and the instruments must be properly cleaned. All the internal and external surfaces must be at least macroscopically free of undesired substances. The cleanliness must be verifiable.

The disinfection procedure must kill vegetative microorganisms and inactivate viruses. The disinfection must be verifiable. All parts which come into contact

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with the material to be processed and the cleaning solution must be disinfected according to the requirements presented above.

### Water Quality

The water should conform to the following quality parameters as nearly as possible and shall be supplied to the process at that quality. Maximum values are:

– Total hardness	5 °d (0.9 mmol/l)
– Chloride	100 mg/l
– Silicic acid/silicate, as SiO <sub>2</sub>	15 mg/l
– Iron	0.05 mg/l
– Manganese	0.05 mg/l
– Copper	0.05 mg/l

Depending on the reprocessing procedure and the agents used, higher values can be tolerated if necessary. Fully demineralised water, which must be of potable water quality with respect to microbiological standards, must be used for the final rinse. Before it is pumped off, this last rinse water must meet the microbiological requirements stated in the (German) Drinking-Water Ordinance. The materials must be free of droplets after the process has terminated.

### Ultrasonic Cleaning

Ultrasound (35 kHz) can be employed in order to assist cleaning the instrument in the joint region and inside the tubes. The effect of ultrasound depends on how the basin is loaded. Exposure should last between 3 to 5 minutes, and the instruments should be placed in the ultrasonic cleaner with their jaws opened.

Which items are suitable for ultrasonic cleaning can be gathered from the manufacturer's data. If the ultrasonic bath is not part of the washer/disinfector, a fume hood over the bath must suck up the vapours produced. During transfer of the items to the washer/disinfector, contamination of the surroundings must be avoided.

### Validation of the Cleaning

The manufacturer of the washer/disinfector must specify a validated standard process for cleaning MIS instruments, including the required chemical components. When the washer/disinfector is installed, the process must be tested on site (qualification and process validation). Changes in the process require revalidation if the new process has not previously been validated.

All processes should be in correspondence with the prospective quality-assurance measures that will shortly be required in hospitals. Deviations from the cleaning programme must be detected by the washer/disinfector, and the unit must automatically cancel any cycle in case of malfunction.

Cleaning and disinfection results must be checked separately after termination of the reprocessing procedure.

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### **Instrument Care and Reassembly**

Because of the complex structure of the instruments and their costly technology, they should be taken apart, reassembled, maintained, and tested for correct function by specially trained personnel only. Inspection and reassembly should be completed before sterilisation. Instruments not cleaned effectively must be processed a second time or removed from use and repaired.

For instrument care, only the joints, sliding surfaces and movable parts are sprayed with an oil suitable for steam sterilisation. Excess oil must be removed with compressed air. Immersion in special baths and spraying on all parts of the instruments is not recommended.

It goes without saying that comprehensive personnel training with regard to the specifics and hazards while processing MIS instruments is a must. Proper and informed handling will also help bring down the high repair bills for these delicate instruments.

### **Sterilisation**

The instruments must be steam-sterilisable and sterilised in assembled condition.

### **Measures in the Operating Room**

Before using an instrument in the operating room, the instrumentation personnel or the surgeon should inspect it one more time. ■

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